




Galaxy II[®]

Self-retaining Surgical Retractor


The Image Guide

Best Surgeon-Focused Medical Device Company 2022
Most Advanced Self-Retaining Surgical Retractor 2022



JUNEMEDICAL
SIMPLIFYING COMPLEXITY

MSM21-7 Rev 12 Last Reviewed MAR 2025. Always read the IFU.
© MOSAIC SURGICAL All rights reserved.


Frames




Product name: Galaxy II Snowman
Product code: JUNE3000



Product name: HandMe Galaxy II
Product code: GX10



Product name: Galaxy II Square
Product code: JUNE3001



Product name: Galaxy II Slider
Product code: JUNEGX-U

All frames sold as:
Box of 10. Frame only, pre-packed & sterile.

Stay Hooks


Product name
Galaxy II 3mm Sharp Hook
Product code
GX40 - Box of 50 individually pouched



Product name
Galaxy II 5mm Sharp Hook
Product codes
JUNE3002s-50 - Box of 50 individually pouched
JUNE3002s-100 - Box of 100 individually pouched




Product name
Galaxy II 5mm Semi-Blunt
Product code
JUNE3002sb-50 - Box of 50 individually pouched



Product name
Galaxy II 7mm Double Sharp Hook
Product code
JUNE3004s-50 - Box of 50 individually pouched



Product name
Galaxy II 5mm Blunt Hook
Product codes
JUNE3002b-50 - Box of 50 individually pouched
JUNE3002b-100 - Box of 100 individually pouched



Product name
Galaxy II 20mm 3-Finger Claw
Product code
JUNE3005b-50 - Box of 50 individually pouched




Product name
Galaxy II 12mm Blunt Hook
Product codes
JUNE3003b-50 - Box of 50 individually pouched
JUNE3003b-100 - Box of 100 individually pouched




Snowman Kits


Pre-packed and sterile



Product name: Galaxy II Snowman Blunt 12mm
Product code: JUNE3000b12
Sold as: Box of 10
Content: 1x Retractor frame
8x 12mm Blunt Hook
Catheter Clip
Light clip



Product name: Galaxy II Snowman Sharp 5mm
Product code: JUNE3000s5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Sharp Hook
Catheter Clip
Light clip




Product name: Galaxy II Snowman Blunt 5mm
Product code: JUNE3000b5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Blunt Hook
Catheter Clip
Light clip




Product name: Galaxy II Snowman Semi-Blunt 5mm
Product code: JUNE3000sb5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Semi-Blunt Hook
Catheter Clip
Light clip

Square Kits


Pre-packed and sterile




Product name: Galaxy II Square Blunt 12mm
Product code: JUNE3001b12
Sold as: Box of 10
Content: 1x Retractor frame
8x 12mm Blunt Hook
Catheter Clip
Light clip



Product name: Galaxy II Square Sharp 5mm
Product code: JUNE3001s5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Sharp Hook
Catheter Clip
Light clip




Product name: Galaxy II Square Blunt 5mm
Product code: JUNE3001b5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Blunt Hook
Catheter Clip
Light clip



Product name: Galaxy II Square Semi-Blunt 5mm
Product code: JUNE3001sb5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Semi-Blunt Hook
Catheter Clip
Light clip



Slider Kits

Pre-packed and sterile




Product name: Galaxy II Slider Blunt 12mm
Product code: JUNE GX-U12b
Sold as: Box of 10
Content: 1x retractor frame
8x 12mm Blunt stays
Catheter clip
Light clip

Product name: Galaxy II Slider Sharp 5mm
Product code: JUNE GX-U5s
Sold as: Box of 10
Content: 1x retractor frame
8x 5mm Sharp stays
Catheter clip
Light clip




Product name: Galaxy II Slider Blunt 5mm
Product code: JUNE GX-U5b
Sold as: Box of 10
Content: 1x retractor frame
8x 5mm Blunt stays
Catheter clip
Light clip

Product name: Galaxy II Slider Semi-Blunt 5mm
Product code: JUNE GX-U5sb
Sold as: Box of 10
Content: 1x retractor frame
8x 5mm Semi-Blunt stays
Catheter clip
Light clip





Galaxy II Male Kit

Product code: JUNE GX-2020
Sold as: Box of 1, pre-packed and sterile
Description: Kit specifically assembled for Male Pelvic Health procedures, for example penile implants.
Content: 1x Snowman Retractor
6x 12mm Blunt Hook
1x 5mm Sharp Hook
2x 3-Finger Claw
2x 7mm Double Sharp Hook
1x Penile Hammock to avoid unnecessary skin puncture
1x Catheter Clip
1x Light clip



How To Galaxy II: Practical Tips




98% of surgeons who try Galaxy II want to switch straight away


100% Surgeon Satisfaction in trial of almost 4000 patients in 12 NHS Trusts

100% CO2 offset every year since launch


Use counter-traction to fasten hooks



Always release hooks from the tissue, never pull on the tail end to release from frame.



History & Recommendations



ENT

FOR MAXIMUM ADJUSTABILITY Order code: JUNE3001sb5
Contents: 1 x Slider, 8 x 5mm SEMI BLUNT stays

Neck

The Galaxy Slider is the most used: Order code: JUNE3001sb5
Contents: 1 x Slider, 8 x 5mm SHARP stays

Hand/arm/elbow

WE USED TO SELL MOSTLY Order code: JUNE3001sb5
Contents: 1 x Square ring, 8 x 5mm SEMI BLUNT stays

But the Galaxy Slider has become a more popular choice.

Dupuytren's or anything in the palm: new HandMe (order code GX30) is specifically designed.

Foot

TRY THIS Order code: JUNE3001sb5
Contents: 1 x Square ring, 8 x 5mm SHARP stays

Plastics

WE RECOMMEND Order code: JUNE3001sb5
Contents: 1 x Slider, 8 x 5mm SHARP stays

Cardio

TRY THIS Order code: JUNE3001sb5
Contents: 1 x Slider, 8 x 5mm BLUNT stays

Hernia

MOST ORDERED Order code: JUNE3000b12(kit) and JUNE3005b(claws)
Contents: 1 x Snowman, 8 x 12mm BLUNT stays (Kit) Single packed additional claws

Male

SPECIFICALLY DESIGNED FOR MALE PROCEDURES Order code: JUNE3001sb5
Content: 6 x 12mm Blunt stays, 1 x 5mm Sharp stay, 2 x 3-finger claw and 2x Double 7mm Sharp Hooks, our unique penile hammock and Snowman shape retractor.

Developed in collaboration with UCLH in London UK

Gyn

MOST POPULAR is Order code: JUNE3000sb5
Contents: 1 x Snowman ring, 8 x 5mm SHARP stays

But....Galaxy Slider with Semiblunt Hook is now available and is tailored specifically for GYN!

Urology


WE RECOMMEND Order code: JUNE3000sb5
Contents: 1 x Snowman ring, 8 x 5mm SEMI BLUNT stays

Colorectal

POPULAR CHOICE is Order code: JUNE3001sb5
Contents: 1 x Square ring, 8 x 5mm BLUNT stays

Paediatrics

Slider kits all the way – so versatile and adjustable to fit. Start with Semiblunt Hooks (sharper or blunter also available) JUNE3001sb5



www.galaxyretractor.com

HandMe® Retractor

a 3-in-1 tool revolutionizing hand surgery



The HandMe Retractor for hand surgery such as Dupuytren's is a new invention specifically designed for improving access and vision for hand surgeons, whilst reducing the clutter and cost of multiple instruments. With its patent pending design, it creates a clean surgical field by covering the fingers and elevating the relevant digits, allowing unprecedented control for the surgeon.

With specially designed digit holders, gentle but very reliable retraction provides both stability and access to the exposed tissue.

The HandMe Retractor replaces the lead hand, the retractor AND the sutures / assistant currently often used in hand surgery: **a 3-in-1 solution!**



HandMe
Galaxy II
surgical retractor

Precision retraction for single or double digits

Reduced tissue trauma

Enhanced visibility and access for surgeons

Ties

There is 80cm of straight, uncut silicone tie included in the HandMe Accessory Pack that can be used to secure additional digits or cut to desired length for holding single digits or light cables to the frame.

The HandMe Retractor consists of a u-shaped half-moon made out of 50gr of medical grade Terlux®. The dimensions are 22cm x 16cm x 3cm. HandMe is designed to fit all patients, both paediatric and adults. The retractor comes sterile and is disposable.



Finger hats

A unique idea to secure digits to the frame without constricting bloodflow, especially important in older patients. Gently curved to firmly hold, with dimpled surface to prevent slipping. Made of Terlux® and silicone. The central part is 2.3cm x 1cm and the tails are 12cm long. There are 2 Finger Hats included in the Accessory Pack, one for single digits and one to hold double digits.

Tissue Hooks

3mm sharp hooks with low profile, made of stainless steel and silicone. There are 6 hooks in the Accessory Pack, each 9cm long.



www.galaxyretractor.com



www.galaxyretractor.com

	Product	Description		Product code
Retractors	HandMe® Galaxy II®			
	HandMe	Frame Only	HandMe® Galaxy II® Retractor Frame Only (Box of 10)	GX10
		Accessories	HandMe® Galaxy II® Accessory Pack 2x Finger Hats, 6x3mm SHARP Stays, 1x Tie (Box of 10)	GX20
		Kit	HandMe® Galaxy II® Retractor Kit (Box of 10) 1x Frame, 2x Finger Hats, 6x3mm SHARP Stays, 1x Tie (Box of 10)	GX30
	Retractor - Galaxy II®			
	Snowman	Frame only	Galaxy II® Snowman Sterile Retractor Frame Only (Box of 10)	JUNE3000
		Kit	Galaxy II® Snowman Sterile Surgical Retractor KIT Frame & 8x5mm SHARP Stays (Box of 10)	JUNE3000s5
			Galaxy II® Snowman Sterile Surgical Retractor KIT Frame & 8x5mm SEMI-BLUNT Stays (Box of 10)	JUNE3000sb5
			Galaxy II® Snowman Sterile Surgical Retractor KIT Frame & 8x5mm BLUNT Stays (Box of 10)	JUNE3000b5
			Galaxy II® Snowman Sterile Surgical Retractor KIT Frame & 8x12mm BLUNT Stays (Box of 10)	JUNE3000b12
	Square	Frame only	Galaxy II® Square Sterile Retractor Frame Only (Box of 10)	JUNE3001
		Kit	Galaxy II® Square Sterile Surgical Retractor KIT Frame & 8x5mm SHARP Stays (Box of 10)	JUNE3001s5
			Galaxy II® Square Sterile Surgical Retractor KIT Frame & 8x5mm SEMI-BLUNT Stays (Box of 10)	JUNE3001sb5
			Galaxy II® Square Sterile Surgical Retractor KIT Frame & 8x5mm BLUNT Stays (Box of 10)	JUNE3001b5
			Galaxy II® Square Sterile Surgical Retractor KIT Frame & 8x12mm BLUNT Stays (Box of 10)	JUNE3001b12
	Slider	Frame only	Galaxy II® Slider Sterile Surgical Retractor Frame Only (Box of 10)	JUNEGX-U
		Kit	Galaxy II® Slider Sterile Surgical Retractor KIT Frame & 8x5mm SHARP Stays (Box of 10)	JUNEGX-U5s
			Galaxy II® Slider Sterile Surgical Retractor KIT Frame & 8x5mm SEMI-BLUNT Stays (Box of 10)	JUNEGX-U5sb
			Galaxy II® Slider Sterile Surgical Retractor KIT Frame & 8x5mm BLUNT Stays (Box of 10)	JUNEGX-U5b
			Galaxy II® Slider Sterile Surgical Retractor KIT Frame & 8x12mm BLUNT Stays (Box of 10)	JUNEGX-U12b
	Male Kit	Kit	Galaxy II® Male KIT 6x12mm Blunt Stays, 1x5mm Sharp Stay, 2x3 Finger Claw and 2x7mm Double Hook Sharp, 1x Penile Hammock and 1x Snowman Retractor Frame	JUNEGX-2020
Stay Hooks	Stay/Hooks - Galaxy II®			
	Sharp 3mm	Individually Packed Hooks	Galaxy II® 3mm SHARP Sterile Surgical Stays (Box of 50)	GX40
	Sharp 5mm		Galaxy II® 5mm SHARP Sterile Surgical Stays (Box of 100)	JUNE3002S-100
			Galaxy II® 5mm SHARP Sterile Surgical Stays (Box of 50)	JUNE3002S-50
	Blunt 5mm		Galaxy II® 5mm BLUNT Sterile Surgical Stays (Box of 100)	JUNE3002B-100
			Galaxy II® 5mm BLUNT Sterile Surgical Stays (Box of 50)	JUNE3002B-50
	Semi-Blunt 5mm		Galaxy II® 5mm SEMI-BLUNT Sterile Surgical Stays (Box of 50)	JUNE3002SB-50
	Blunt 12mm		Galaxy II® 12mm BLUNT Sterile Surgical Stays (Box of 100)	JUNE3003B-100
			Galaxy II® 12mm BLUNT Sterile Surgical Stays (Box of 50)	JUNE3003B-50
	Double Sharp 7 mm		Galaxy II® 7mm Double Hook SHARP Sterile Surgical Stays (Box of 50)	JUNE3004S-50
	3 Finger Claw		Galaxy II® 3 Finger CLAW Sterile Surgical Stays (Box of 50)	JUNE3005B-50

