



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2021



Galaxy II[®]

Self-retaining Surgical Retractor

The Image Guide

Best Surgeon-Focused Medical Device Company 2022
Most Advanced Self-Retaining Surgical Retractor 2022: Galaxy II[®]



JUNEMEDICAL
SIMPLIFYING COMPLEXITY

MSM21-7 Rev 9 Last Reviewed JUN 2022. Always read the IFU. © MOSAIC SURGICAL All rights reserved.

Frames



Product name: Galaxy II Snowman
Product code: JUNE3000
Sold as: Box of 10
Description: Frame only
Pre-packed and sterile
Size: Inner Length 32cm (13 inches)
Large Inner Diameter 17.78cm (7 inches)
Small Inner Diameter 11.2cm (4.4 inches)



Product name: Galaxy II Square
Product code: JUNE3001
Sold as: Box of 10
Description: Frame only
Pre-packed and sterile
Size: Inner Diameter 14x14cm (5.5x 5.5 inches)



Product name: Galaxy II Slider
Product code: JUNEGX-U
Sold as: Box of 10
Description: Frame only
Pre-packed and sterile
Size: Inner length 13cm (5.12 inches)
Inner diameter when closed 11.3 cm (4.45 inches)
Inner diameter when fully open 16cm (6.30 inches)

Stay Hooks


Product name: Galaxy II Blunt Small Hook

Size: 5mm

Product code: JUNE3002b

Sold as: Box of 50 individually pouched (JUNE3002b-50)

Box of 100 individually pouched (JUNE3002b-100)



Product name: Galaxy II Sharp Small 5mm Hook

Size: 5mm

Product code: JUNE3002s

Sold as: Box of 50 individually pouched (JUNE3002s-50)

Box of 100 individually pouched (JUNE3002s-100)




Product name: Galaxy II Blunt Large Hook

Size: 12mm

Product code: JUNE3003b

Sold as: Box of 50 individually pouched (JUNE3003b-50)

Box of 100 individually pouched (JUNE3003b-100)



Product name: Galaxy II Double Sharp Hook

Size: 7mm

Product code: JUNE3004s-50

Sold as: Box of 50 individually pouched




Product name: Galaxy II 3-Finger Claw

Size: 20mm

Product code: JUNE3005b-50

Sold as: Box of 50 individually pouched



Product name: Galaxy II Semi-Blunt Small

Size: 5mm

Product code: JUNE3002sb-50

Sold as: Box of 50 individually pouched




 **GALAXY II**

www.galaxyretractor.com

Snowman Kits

Pre-packed and sterile



Product name: Galaxy II Snowman Blunt 12mm

Product code: JUNE3000b12

Sold as: Box of 10

Content: 1x Retractor frame

8x 12mm Blunt Hook

Catheter Clip

Galaxy LUXC clip for light attachment

Product name: Galaxy II Snowman Sharp 5mm

Product code: JUNE3000s5


Sold as: Box of 10

Content: 1x Retractor frame

8x 5mm Sharp Hook

Catheter Clip

Galaxy LUXC clip for light attachment



Product name: Galaxy II Snowman Blunt 5mm

Product code: JUNE3000b5

Sold as: Box of 10

Content: 1x Retractor frame

8x 5mm Blunt Hook

Catheter Clip

Galaxy LUXC clip for light attachment

Product name: Galaxy II Snowman Semi-Blunt 5mm

Product code: JUNE3000sb5



Sold as: Box of 10

Content: 1x Retractor frame

8x 5mm Semi-Blunt Hook

Catheter Clip

Galaxy LUXC clip for light attachment



Need light?

Add -LUXC for light accessory

e.g. JUNE3000B5-LUXC




GALAXY II

www.galaxyretractor.com



Square Kits

Pre-packed and sterile





Product name: Galaxy II Square Blunt 12mm
Product code: JUNE3001b12
Sold as: Box of 10
Content: 1x Retractor frame
8x 12mm Blunt Hook
Catheter Clip
Galaxy LUXC clip for light attachment

Product name: Galaxy II Square Sharp 5mm
Product code: JUNE3001s5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Sharp Hook
Catheter Clip
Galaxy LUXC clip for light attachment



Product name: Galaxy II Square Blunt 5mm
Product code: JUNE3001b5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Blunt Hook
Catheter Clip
Galaxy LUXC clip for light attachment

Product name: Galaxy II Square Semi-Blunt 5mm
Product code: JUNE3001sb5
Sold as: Box of 10
Content: 1x Retractor frame
8x 5mm Semi-Blunt Hook
Catheter Clip
Galaxy LUXC clip for light attachment



Need light?


Add -LUXC for light accessory
e.g. JUNE3001B5-LUXC

 **GALAXY II**

www.galaxyretractor.com



Slider Kits

Pre-packed and sterile




Product name: Galaxy II Slider Blunt 12mm
Product code: JUNEGX-U12b
Sold as: Box of 10
Content: 1x retractor frame
8x 12mm Blunt stays
Catheter clip
Galaxy LUXC clip for light attachment

Product name: Galaxy II Slider Sharp 5mm
Product code: JUNEGX-U5s
Sold as: Box of 10
Content: 1x retractor frame
8x 5mm Sharp stays
Catheter clip
Galaxy LUXC clip for light attachment



Product name: Galaxy II Slider Blunt 5mm
Product code: JUNEGX-U5b
Sold as: Box of 10
Content: 1x retractor frame
8x 5mm Blunt stays
Catheter clip
Galaxy LUXC clip for light attachment

Product name: Galaxy II Slider Semi-Blunt 5mm
Product code: JUNEGX-U5sb
Sold as: Box of 10
Content: 1x retractor frame
8x 5mm Semi-Blunt stays
Catheter clip
Galaxy LUXC clip for light attachment



Need light?

Add -LUXC for light accessory
e.g. JUNE3001B5-LUXC

Galaxy II Male Kit

Product code: JUNEGX-2020

Sold as: Box of 1, pre-packed and sterile

Description: Kit specifically assembled for Male Pelvic Health procedures, for example penile implants.

Content: 1x Snowman Retractor

6x 12mm Blunt Hook



1x 5mm Sharp Hook

2x 3-Finger Claw

2x 7mm Double Sharp Hook

1x Penile Hammock to avoid unnecessary skin puncture

Catheter Clip and LUXC Clip for light attachment



Need light?

Add -LUXC for light accessory
e.g. JUNEGX-2020-LUXC

 **GALAXY II**

www.galaxyretractor.com

Product Code Overview

Product		Description	Product Code	
Retractor Kits	Snowman	Frame Only	GALAXY II™ Snowman Sterile Retractor Frame Only (Box of 10)	JUNE3000
		Kit	GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SHARP (Box of 10)	JUNE3000s5
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT (Box of 10)	JUNE3000sb5
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm BLUNT (Box of 10)	JUNE3000b5
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 12mm BLUNT (Box of 10)	JUNE3000b12
	Square	Frame Only	GALAXY II™ Square Sterile Retractor Frame Only (Box of 10)	JUNE3001
		Kit	GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SHARP (Box of 10)	JUNE3001s5
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT (Box of 10)	JUNE3001sb5
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm BLUNT (Box of 10)	JUNE3001b5
			GALAXY II™ Square Sterile Surgical Retractor KIT - 12mm BLUNT (Box of 10)	JUNE3001b12
	Slider	Frame Only	GALAXY II™ Slider Sterile Surgical Retractor Frame Only (Box of 10)	JUNEGX-U
		Kit	GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SHARP (Box of 10)	JUNEGX-U5s
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT (Box of 10)	JUNEGX-U5sb
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm BLUNT (Box of 10)	JUNEGX-U5b
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 12mm BLUNT (Box of 10)	JUNEGX-U12b
	Male Kit	Kit	6x 12mm Blunt Stays, 1x 5mm Sharp Stay, 2x 3 Finger Claw and 2x 7mm Double Hook Sharp, 1x Penile Hammock and 1x Snowman Retractor Frame	JUNEGX-2020
Stay Hooks	Sharp 5mm	Individually Packed Hooks	GALAXY II™ 5mm SHARP Sterile Surgical Stays (Box of 50)	JUNE3002s-50
	Semi-Blunt 5mm		GALAXY II™ 5mm SHARP Sterile Surgical Stays (Box of 100)	JUNE3002s-100
	Blunt 5mm		GALAXY II™ 5mm SEMI-BLUNT Sterile Surgical Stays (Box of 50)	JUNE3002sb-50
			GALAXY II™ 5mm BLUNT Sterile Surgical Stays (Box of 50)	JUNE 3002b-50
	Blunt 12mm		GALAXY II™ 5mm BLUNT Sterile Surgical Stays (Box of 100)	JUNE 3002b-100
			GALAXY II™ 12mm BLUNT Sterile Surgical Stays (Box of 50)	JUNE 3003b-50
	Double Sharp 7 mm		GALAXY II™ 12mm BLUNT Sterile Surgical Stays (Box of 100)	JUNE 3003b-100
3 Finger Claw	GALAXY II™ 7mm Double Hook SHARP Sterile Surgical Stays (Box of 50)	JUNE3004s-50		
		GALAXY II™ 20mm 3 Finger Claw BLUNT Sterile Surgical Stays (Box of 50)	JUNE3005b-50	
Retractor Kits With Light	Snowman LUX Connect	Frame Only	GALAXY II™ Snowman Sterile Retractor and Sterile Surgical Light (Box of 10)	JUNE3000-LUXC
		Kit	GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SHARP Stays and Sterile Surgical Light (Box of 10)	JUNE3000s5-LUXC
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3000sb5-LUXC
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3000b5-LUXC
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 12mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3000b12-LUXC
	Square LUX Connect	Frame Only	GALAXY II™ Square Surgical Retractor Frame and Sterile Surgical Light (Box of 10)	JUNE3001-LUXC
		Kit	GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SHARP Stays and Sterile Surgical Light (Box of 10)	JUNE3001s5-LUXC
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3001sb5-LUXC
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3001b5-LUXC
			GALAXY II™ Square Sterile Surgical Retractor KIT - 12mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3001b12-LUXC
	Slider LUX Connect	Frame Only	GALAXY II™ Slider Sterile Surgical Retractor Frame and Sterile Surgical Light (Box of 10)	JUNEGX-U-LUXC
		Kit	GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SHARP Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U5s-LUXC
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U5sb-LUXC
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U5b-LUXC
GALAXY II™ Slider Sterile Surgical Retractor KIT - 12mm BLUNT Stays and Sterile Surgical Light (Box of 10)			JUNEGX-U12b-LUXC	
Male Kit LUX Connect	Kit	GALAXY II™ Male KIT 6x 12mm Blunt Stays, 1x 5mm Sharp Stay, 1x 3 Finger Claw and 2x 7mm Double Sharp Hooks, 1x Penile Hammock, 1x Snowman Retractor Frame and 1x Sterile Surgical Light	JUNEGX-2020-LUXC	
Light	JUNELUX Connect	Single Light	LUX Connect Sterile Surgical Light	JUNELUXC1

GALAXY II

The 8 points of retraction



How do you prefer your surgical site to look?

Order from: order@junemedical.com

The problems

After observing 100's of hours of surgeries, we started to notice some very common, recurring problems facing surgeons and medical staff.



Loose parts

Loose parts are always a major issue in any surgical environment. They heighten the risk and increase unnecessary work for staff.



Poor construction

Old screw mechanisms led to retractors not tightening enough, while over-tightening caused screws to break and fall off.



Bad quality

Many of the models on today's current market are very poorly constructed and do not have the required certification.



Outdated design

Although revolutionary 15 years ago, in today's medical environment they are outdated, clunky and invasive.



Too expensive

Even though these models have not changed in over a decade, they remain very expensive. Design, materials and manufacturing processes have improved dramatically, and associated costs can be reduced.



The Galaxy II provides 8 points of retraction, while I only have two hands!

Kriston Bertelsen - Surgical Tech.
Four Peaks Surgery Center AZ,
United States





The solutions

When you introduce the Galaxy II, those identified common problems become things of the past. The Galaxy II provides clear visualization, optimal exposure and access to the operative area during surgical procedures. The first generation Galaxy, invented in 2016, assisted many surgical procedures, but we knew there was room for improvement.

So we attended several hours of surgery, asked lots of questions and collated a multitude of feedback.

The result is the award winning Galaxy II.



Cam locks

Intuitive to use, the surgeon inspired CAM locking system enables users to swiftly adapt the retractor to the needs of the procedure with one single hand.



Lightweight

The cored out frame is made of a strong material, and is still 10% lighter than it's predecessor.



Precise design

The slots on the ring easily secure the stays in place, but readily release when rapid reposition is required. The low profile allows superior access to the operative sites.



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2021



JUNEMEDICAL

SIMPLIFYING COMPLEXITY



www.junemedical.com



[junemedical](https://twitter.com/junemedical)



contact@junemedical.com



[june-medical](https://www.linkedin.com/company/june-medical)

Galaxy II is distributed by JUNE Medical

JUNE Medical Ltd, 9 York Way, High Wycombe, Bucks, HP12 3PY, UK
JUNE Medical USA Inc. 257 Old Churchmans Road, New Castle, DE 19720 USA
JUNE Medical Sverige AB, Alsbergsvagen 10, 51156 Kinna, Sweden



JUNE Medical Sverige AB
Alsbergsvagen 10
51156 Kinna
Sweden



Galaxy II is manufactured by Mosaic Surgical Limited

9 York Way, High Wycombe, Bucks, HP12 3PY, UK
contact@mosaicsurgical.com