



# GalaxyII<sup>®</sup>

Self-retaining Surgical Retractor

Product Brochure

Best Surgeon-Focused Medical Device Company 2022  
Most Advanced Self-Retaining Surgical Retractor 2022: Galaxy II<sup>®</sup>

# The problems

After observing 100's of hours of surgeries, we started to notice some very common, recurring problems facing surgeons and medical staff.



## Loose parts

Loose parts are always a major issue in any surgical environment. They heighten the risk and increase unnecessary work for staff.



## Poor construction

Old screw mechanisms led to retractors not tightening enough, while over-tightening caused screws to break and fall off.



## Bad quality

Many of the models on today's current market are very poorly constructed and do not have the required certification.



## Outdated design

Although revolutionary 15 years ago, in today's medical environment they are outdated, clunky and invasive.



## Too expensive

Even though these models have not changed in over a decade, they remain very expensive. Design, materials and manufacturing processes have improved dramatically, and associated costs can be reduced.



The Galaxy II provides 8 points of retraction, while I only have two hands!

**Kriston Bertelsen** - Surgical Tech.  
Four Peaks Surgery Center AZ,  
United States





# The solutions

When you introduce the Galaxy II, those identified common problems become things of the past. The Galaxy II provides clear visualization, optimal exposure and access to the operative area during surgical procedures. The first generation Galaxy, invented in 2016, assisted many surgical procedures, but we knew there was room for improvement.

So we attended several hours of surgery, asked lots of questions and collated a multitude of feedback.

The result is the award winning Galaxy II.



## Cam locks

Intuitive to use, the surgeon inspired CAM locking system enables users to swiftly adapt the retractor to the needs of the procedure with one single hand.



## Lightweight

The cored out frame is made of a strong material, and is still 10% lighter than it's predecessor.



## Precise design

The slots on the ring easily secure the stays in place, but readily release when rapid reposition is required. The low profile allows superior access to the operative sites.



## The additional benefits

A self-retaining, single use retractor system providing clear visualisation, optimal exposure and access to the operative area.



### Made in Britain

We manufacture Galaxy II in the UK to ensure great quality in every aspect. We bring only the best medical devices to our customers.



### Swedish design

Sleek, light, clean and innovative design - optimised for safety, efficiency and manufacturing excellence.



### The popular choice


**98% of surgeons** who trial respond that they will CONTINUE to use Galaxy II. Data on file.



Department for  
International Trade



Galaxy II is a very useful device, and I use it for all my cases where I need effective retraction and good access. Good design improvements over the older retractors.”




...really light weight which makes a big difference. The locks are good: they really make it stay tight – I could never get the old screws to really do that!”




# Frames and accessories


Combine any way you want for any type of surgery; ENT, General, Orthopaedic, Urology, Gynaecology, Breast, Plastic, to name a few.



## 3x Frames



## 6x Stays




# Product Order Codes

	Product	Description	Product Code
Stay Hooks	Sharp 5mm	GALAXY II™ 5mm SHARP Sterile Surgical Stays (Box of 50)	JUNE3002s-50
		GALAXY II™ 5mm SHARP Sterile Surgical Stays (Box of 100)	JUNE3002s-100
	Semi-Blunt 5mm	GALAXY II™ 5mm SEMI-BLUNT Sterile Surgical Stays (Box of 50)	JUNE3002sb-50
		GALAXY II™ 5mm BLUNT Sterile Surgical Stays (Box of 50)	JUNE 3002b-50
	Blunt 5mm	GALAXY II™ 5mm BLUNT Sterile Surgical Stays (Box of 100)	JUNE 3002b-100
		GALAXY II™ 12mm BLUNT Sterile Surgical Stays (Box of 50)	JUNE 3003b-50
	Blunt 12mm	GALAXY II™ 12mm BLUNT Sterile Surgical Stays (Box of 100)	JUNE 3003b-100
	Double Sharp 7 mm	GALAXY II™ 7mm Double Hook SHARP Sterile Surgical Stays (Box of 50)	JUNE3004s-50
	3 Finger Claw	GALAXY II™ 20mm 3 Finger Claw BLUNT Sterile Surgical Stays (Box of 50)	JUNE3005b-50


# Pre-packed kits

Simple, cost-effective, space-saving and environmentally friendly with 10% less plastic. Our pre-packed and sterile single-use kits are the clear winners across almost all specialities.



## Square Kit

Square frame with blunt, sharp and semi-blunt hooks.



## Snowman Kit

Snowman frame with blunt, sharp and semi-blunt hooks.



## Slider Kit

Slider frame with blunt, sharp and semi-blunt hooks.

# JUNEGX-2020 Male Kit

Specifically designed to meet the needs of male urology and andrology procedures. Designed with inspiration from UCLH in London, UK, featuring innovative solutions to reduce tissue trauma and risk of infection using a novel penile hammock.



Light weight



Pre-packed for simplicity



Strong and reliable



Greater visual area




Innovative cam locks



Globally approved

## Product Order Codes

Product	Description	Product Code		
Retractor - Galaxy II				
Retractor Kits	Snowman	Frame Only GALAXY II™ Snowman Sterile Retractor Frame Only (Box of 10)	JUNE3000	
		Kit	GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SHARP (Box of 10)	JUNE3000s5
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT (Box of 10)	JUNE3000sb5
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm BLUNT (Box of 10)	JUNE3000b5
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 12mm BLUNT (Box of 10)	JUNE3000b12
	Square	Frame Only GALAXY II™ Square Sterile Retractor Frame Only (Box of 10)	JUNE3001	
		Kit	GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SHARP (Box of 10)	JUNE3001s5
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT (Box of 10)	JUNE3001sb5
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm BLUNT (Box of 10)	JUNE3001b5
			GALAXY II™ Square Sterile Surgical Retractor KIT - 12mm BLUNT (Box of 10)	JUNE3001b12
	Slider	Frame Only GALAXY II™ Slider Sterile Surgical Retractor Frame Only (Box of 10)	JUNEGX-U	
		Kit	GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SHARP (Box of 10)	JUNEGX-U5s
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT (Box of 10)	JUNEGX-U5sb
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm BLUNT (Box of 10)	JUNEGX-U5b
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 12mm BLUNT (Box of 10)	JUNEGX-U12b
	Male Kit	Kit GALAXY II™ Male KIT 6x 12mm Blunt Stays, 1x 5mm Sharp Stay, 2x 3 Finger Claw and 2x 7mm Double Hook Sharp, 1x Penile Hammock and 1x Snowman Retractor Frame	JUNEGX-2020	



# GALAXY II LUX CONNECT

## The world's 1st ring retractor with a light

Single use, continuous light powered by a standard light source.

The LUX Connect is a minimally invasive, dual-lead design that provides light coverage at two points, removing shadows and fully exposing the surgical site to optimise patient safety.


The LUX Connect attaches to the Galaxy II frame via custom-made LUX clips and can be positioned in any of the frame slots.

	Product	Description	Product Code	
Retractor Kits With Light	Snowman LUX Connect	Frame Only	GALAXY II™ Snowman Sterile Retractor and Sterile Surgical Light (Box of 10) JUNE3000-LUXC	
		Kit	GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SHARP Stays and Sterile Surgical Light (Box of 10)	JUNE3000s5-LUXC
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3000sb5-LUXC
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 5mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3000b5-LUXC
			GALAXY II™ Snowman Sterile Surgical Retractor KIT - 12mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3000b12-LUXC
	Square LUX Connect	Frame Only	GALAXY II™ Square Surgical Retractor Frame and Sterile Surgical Light (Box of 10) JUNE3001-LUXC	
		Kit	GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SHARP Stays and Sterile Surgical Light (Box of 10)	JUNE3001s5-LUXC
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3001sb5-LUXC
			GALAXY II™ Square Sterile Surgical Retractor KIT - 5mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3001b5-LUXC
			GALAXY II™ Square Sterile Surgical Retractor KIT - 12mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNE3001b12-LUXC
	Slider LUX Connect	Frame Only	GALAXY II™ Slider Sterile Surgical Retractor Frame and Sterile Surgical Light (Box of 10) JUNEGX-U-LUXC	
		Kit	GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SHARP Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U5s-LUXC
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm SEMI-BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U5sb-LUXC
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 5mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U5b-LUXC
			GALAXY II™ Slider Sterile Surgical Retractor KIT - 12mm BLUNT Stays and Sterile Surgical Light (Box of 10)	JUNEGX-U12b-LUXC
	Male Kit LUX Connect	Kit	GALAXY II™ Male KIT 6x 12mm Blunt Stays, 1x 5mm Sharp Stay, 1x 3 Finger Claw and 2x 7mm Double Sharp Hooks, 1x Penile Hammock, 1x Snowman Retractor Frame and 1x Sterile Surgical Light JUNEGX-2020-LUXC	
Light	JUNELUX Connect	Single Light	LUX Connect Sterile Surgical Light JUNELUXC1	



# GALAXY II

## The 8 points of retraction



How do you prefer your surgical site to look?

Order from: [order@junemedical.com](mailto:order@junemedical.com)



THE QUEEN'S AWARDS  
FOR ENTERPRISE:  
INNOVATION  
2021



# JUNEMEDICAL

SIMPLIFYING COMPLEXITY



[www.junemedical.com](http://www.junemedical.com)



[junemedical](https://twitter.com/junemedical)



[contact@junemedical.com](mailto:contact@junemedical.com)



[june-medical](https://www.linkedin.com/company/june-medical)

**Galaxy II is distributed by JUNE Medical**

JUNE Medical Ltd, 9 York Way, High Wycombe, Bucks, HP12 3PY, UK  
JUNE Medical USA Inc. 257 Old Churchmans Road, New Castle, DE 19720 USA  
JUNE Medical Sverige AB, Alsbergsvagen 10, 51156 Kinna, Sweden



JUNE Medical Sverige AB  
Alsbergsvagen 10  
51156 Kinna  
Sweden



**Galaxy II is manufactured by Mosaic Surgical Limited**

9 York Way, High Wycombe, Bucks, HP12 3PY, UK  
[contact@mosaicsurgical.com](mailto:contact@mosaicsurgical.com)