

# Haemoband-PLUS

Haemorrhoid Multi-Band Ligator



Description	Dimension	Material
External Nozzle End	Ø 15mm	Terlux 2812 HD Polymer
Internal Nozzle End	Ø 10mm	Terlux 2812 HD Polymer
Latex-Free Bands	Internal - Ø 2mm External - Ø 6mm	Polysisoprene 266F10GAS5



Automatic Reloading Technology

Single-Use

## Clinically Efficient

The multi-action handle requires only one hand to control suction, release and reload bands.

## Simple & Fast

Pre-loaded with 4 bands on see-through nozzle. The procedure is reduced from 20 minutes to 2-4 mins.

## Superior Banding

The individually moulded bands roll off very easily, eliminating twisting and ensuring a reliable release and secure fit onto tissue.

## Ready To Use

From an airtight, clean bag.

## Cost Effective

The most cost-effective multi-band ligators available.



haemobandsurgical.com