

Universal  
Applicator

made in  
Germany

# HaemoCer<sup>TM</sup>

## Universal Applicator

For all minimally invasive surgeries

Fast - Effective - Safe

HaemoCer<sup>TM</sup> Universal Applicator is designed to deliver the haemostatic HaemoCer<sup>TM</sup> powder in laparoscopic, ENT, spinal and other MIS procedures. It guarantees precise application of APH particles to the bleeding site where this is mandatory. The delivery instrument can be easily attached to the HaemoCer<sup>TM</sup> or HaemoCer<sup>TM</sup> PLUS bellows dispenser and enables hemostat delivery under direction to the wound site for the control of capillary, venous and arteriole bleeding.



440 mm applicator



High flexibility



Integrated closure clip




Innovative applicator tip

BioCer



# HaemoCer™

## Universal Applicator



- Application in minimally invasive surgical procedures
- Longer applicator for deeper sourced bleeding
- Innovative applicator tip features a design which minimizes risks of blockage in case of contact with moisture
- Integrated shut-off clamp prevents the powder from uncontrolled discharge
- Highly flexible enabling precise tip placement and offering enhanced handling and maximizing wound coverage capability
- Push-on coupling sheath offering simple and firm connection to HaemoCer™ or HaemoCer™ PLUS



### Product Overview HaemoCer™ Universal Applicator

HFZ101

440 mm

5 pcs

# BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0  
Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de

